648459

ko0225

EN Massager gun - User manual

U head / for vertebrae

The reasonable U-shaped design perfectly delivers force path between two pressure points, which is even and powerful. Suitable for stimulating massage such as neck and Spine.

Small round head/ for joints

Small round head brings a powerful yet comfortable pressure for your muscle groups. Suitable for muscle relaxation in all muscle area.

Ball head / for large muscle groups

Ball head design produces a powerful output. Soft material brings an even pressure to your muscle. It is suitable for any large muscle groups over your body.

Small flat head / for all parts of the body

Small flat head is small but full of impact and brings direct pressure to the depths. Great for deep tissues such as palm, sole of foot and joints.

Usage instructions

The deep muscle fascia massage gun is a portable and rechargeable vibration device with large torque and different speed adjustment levels. It can effectively reduce the pain caused by the accumulation of lactic acid after intense exercise. Enhance your well-being by shaping the muscle and relieving muscle discomfort.

Instructions

1. Power on / off button (long press to start / long press to turn off).
2. Press the key gear to change.
3. Detachable and replaceable vibrator. When disassembling the vibrator, just pull it out with force.
4. Please select 5V charging equipment for charging.
5. For the first time, please charge for 3 hours. Please connect the data cable to the 5V adapter for charging.
6. Button type: when charging, 6 lights will flash, when full, 6 lights will be on, and when working, the lights will indicate the gear position.
7. Touch type: when charging, the nixie tube displays the electric quantity, while working, it displays the gear.

Suggested application

Attach the head to the device, Safely power on the device before contact with your body. Never use the device when charging to avoid risks. Move the device along the entire length of the muscle at a moderate pace for 30s to activate the muscle. To relax the muscle, move the device slowly along the entire length of your muscle for 2 minutes maximum each time. For best results, use light to moderate pressure when applying the device. Repeat this regime for different areas of muscles. Please note that this is only a suggested application for the use of massage gun. Please be aware that bruising may occur. Result may vary from person to person.

Safety instructions

Danger: In order to reduce the risk of electric shock, fire and personal injury, this product must be used in accordance with the following instructions.

- For adults only
- It can only be used on dry, clean skin surfaces of the body, and cannot be used through clothes. Press and move it gently on the skin. The duration on each part is about 60 seconds.

- Use this device only on soft tissues of the body without causing pain or discomfort. Do not use it on head or any hard or boned part of the body.
- Only use the vibrating head that can produce the best results.
- Bruise may be caused whether under moderate pressure or normal pressure. Please always check the area that receives physiotherapy.
- Stop using it immediately once there is any sign of pain or discomfort.
- Keep fingers, hair or other parts of the body away from the shaft and back of the vibrating head, because they might be pinched.
- Do not place any object in the vent of the device.
- Do not immerse it in water or allow water to enter the vent of the device.
- Do not drop the device or use it improperly.
- Charge the device charger with only 8.4 VDC.
- Carefully check the device and battery before each use.
- Do not change or modify the device in any way.
- Never operate or charge the device in an unattended way.

Disclaimer

PLEASE DO NOT USE THE MASSAGE GUN OR ANY VIBRATION DEVICE WITHOUT FIRST OBTAINING APPROVAL FROM YOUR DOCTOR IF ANY OF THE FOLLOWING APPLY:

Pregnancy, diabetes with complications such as neuropathy or retinal damage. Wear of pacemakers, recent surgery, epilepsy or migraines, herniated disks, spondylolysis, recent joint replacement, IUD, meralgia or any concerns about your physical health. Frail individuals and children should be accompanied by a healthy adult when using any vibration device. These contra indications do not mean that you are not able to use a vibration or exercise device, but we advise you to consult a doctor first. Ongoing research is being done to the effects of vibration exercise for specific medical disorders. This will likely result in the shortening of the list of contra indications as shown above. Practical experience has shown the integration of vibration exercise into a treatment plan is advisory in a number of cases. This must be done on the advice of and in accompaniment with a doctor, specialist or physiotherapist.

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