EN Body Fat Scale ZT5106B / 601349

Manual

Technical data

Capacity: 3–180 kg / 6.6–396 lb

Display: 4-digit LCD

Battery: 3 V/2*AAA

Operating temperature: 10–35 °C

Unit: kg/lb/st:lb

Fat parameter display mode: fat percentage, moisture percentage, muscle content, skeletal proportions, calorie recommendation

Auto: weight of the foot

Switch off: after 15 seconds

Low battery indication: "Lo"

Error indication: "Err"

Operation

ACTIVATION FROM THE STANDBY MODE

1) Only weight 2) Set user parameters 3) Measurement of body fat

Notice: The switching operation can only be carried out in the standby mode. Choose the unit "KG", "LB" or "ST" before the operation.

1) Only weight:

Place the scale on a solid level floor. Step on the scale steadily. The figure showing your weight will flash and be locked. The figure will remain on the LCD for 15 seconds, then the scale switches off.

2) Set user parameters:

1. User ID setting: In the standby state, press the "SET" key. The display will show "P0" flashing. Press the arrow key "DOWN" or "UP" to select the one of the 10 user IDs (P0–P9).
2. Height setting: Press the "SET" key to confirm the gender and enter the height setting mode. Press the arrow key "DOWN" or "UP" to set the height (100–220 cm). You can press the "SET" key to confirm the height.
3. Age setting: Press the "SET" key to confirm the height and enter the age setting mode. Press the arrow key "DOWN" or "UP" to set the age (18–100). You can press the "SET" key to confirm the age.
4. Gender setting: Press the "SET" key to confirm the user ID and enter the gender setting mode. Press the arrow key "DOWN" or "UP" to choose the male or female symbol. The LCD will show "0.00". All settings are displayed for 2 seconds.

3) Measurement of body fat:

When the LCD shows "0.0", step barefoot on the scale, keep the sole of your feet on the two electrode plates in full contact. Stand steadily in the middle of scale and start the weighing. After weighing, the figure on the LCD showing your weight will stop flashing and be locked. The weight will be displayed for 3 seconds. Then the body fat percentage, body water content, somatic muscle content, skeletal proportions and calorie recommendation will each be displayed for 3 seconds. This is repeated twice, and then the scale will switch off automatically.

Attention

Avoid loading something on the scale if you leave it unused so as to save battery power.

Keep your feet in the middle of the scale to avoid slipping.

Keep your body steady on the scale so that your weight is measured correctly and smoothly.

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